DES KROKODILS NÄCHSTE GENERATION
THE CROCODILE'S NEXT GENERATION

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Since its introduction at the 1995 IDS in Cologne, Germany, Vision Classic has established itself as an easy to use ceramic with consistent shade accuracy. As a clinically proven ceramic, Vision Classic is now used in over 36 countries.

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*Compatible with alloys having a CTE of 13.8 up to 15.2µm/mK (25 - 600°C)**

- o Precious alloys
- o Semi-precious
- o Non-precious
- o Non greening in combination with alloys containing silver

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
Wohlwend's unique processing yields ideal grain size and pigment distribution. Benefits of this unique process include:

- o Consistently accurate shade reproduction - V shades A1 - D4
- o Easy handling
- o No slumping during build up
- o Minimal
- o ***Less* shrinkage during firing**

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The vital shade reproduction and natural fluorescence of Vision porcelains correspond to natural teeth. This replication of nature can be obtained by either using a basic two layer buildup or by incorporating additional dentine and incisal modifiers into the restoration.

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*Metal processing prior
to porcelain application*

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
Process alloys according to manufacturers instructions (Delete "and avoid sharp edges"), Unless specified, the alloy should be sandblasted with min. 110µm AL²O² at 2 Bar.

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
In most cases, especially on NP, the oxides are to be removed using min. 110µm aluminum oxide. (Refer to manufacturers manual!) A clean and grease-free surface must be present.

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*Apply opaque with a 75-80% semi-
covering layer.*

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*Refer to firing chart.
Shiny appearance after firing.*

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*Powder / Liquid: Mix to creamy consistan-
cy Apply using brush or glass instument.
Firing temperature 970°C*

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*Vis Opaque - Paste Apply using brush or
glass instument. Firing temperature
970°C Can be reduced to degas tempera-
ture of alloys. Prolong hold time at final
temp. 2-3 min.*

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

To achieve internal effects, characterizations can be made by using opaque-modifiers "IO"

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
To stabilize the shade in thin areas, apply a layer of opacious dentine.

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Incisal / Transpa (Optional. To optimize light transmission)

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To achieve opalescence in the incisal area, TR Opal can be mixed into all transpa materials. Even on a simple build up, natural light transmission and reflection can be achieved.

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
Mix glaze powder with liquid and apply in a thin layer. For individual characterisation use stains.

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Depending on desired level of translucency, apply TR Clear or Super Clear.

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Glaze appearance can be controlled by variation of hold time on final temperature.

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By using our range of further materials and simple steps of layering , a natural looking reconstruction can be created. These can be used straight or mixed into other materials.

Erweitert / Modified	Farbe / Colour	A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
	Schultermassen / Shoulder material	SP Neutral + BA					SP Neutral + BB				SP Neutral + BC				SP Neutral + BD		
Opaker zur Chromasteigerung <i>Opaque to intensify chroma</i>	OBA						O BB				O BC				OBD		
		10%	20%	30%	35%	40%	10%	20%	30%	40%	10%	20%	30%	40%	20%	30%	40%
Intensivopaker zum charakterisieren <i>Intensiv opaques for characterization</i>	I O W GB OR OK BR V																
Dentin zur Chromasteigerung <i>Body to intensify chroma</i>	DBA						D BB				D BC				DDB		
		10%	20%	30%	35%	40%	10%	20%	30%	40%	10%	20%	30%	40%	20%	30%	40%
Intensivdentine zum charakterisieren <i>Intensiv dentine for characterization</i>	I D EB ST HG BR MG BL NEON																
Intensivschneide zum individualisieren <i>Intensiv incisal for individualization</i>	I S 7 8 9 10																

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

Reduce dentin slightly more than usual

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*Incisal S58 and transpa TR Opal are alternately
layered to achieve the shade but also to have
natural opaleszenz.*

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*TR Opal can be mixed with all transpas,
depending on which level of translucency is
to be achieved in correction*

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
*To raise the chroma mix 30% ID BA toward
dentine A3.*

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*Mix the highly fluorescent ID Neon toward
the intensive incal IS 7, to achieve higher
value and apply in a thin layer.*

*Mix ID Neon and stain corresponding to mamelon
appearance with stain liquid.
wash onto wet surface.
(What you see is what you get!)*

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Firing chart PFM

Recommended Firing parameters!

Temperatures can vary due to furnace performance.

	Start temperature	Dry time	Heat rate	Final temperatur	Hold	Vacuum	Appearanc
Degassing	Follow alloy manufacturers manual						
Opaque 1	1112 °F	2 Min.	176 °F/Min.	1778 °F	1 Min	Yes	Shiny
Opaque 2	1112 °F	3 Min	176 °F/Min.	1742 °F	1 Min	Yes	Shiny
Shoulder 1 / 2	842 °F	3 Min	113 °F/Min.	1724 °F	1 Min	Yes	Shiny
Dentine 1	842 °F	6 Min	113 °F/Min	1688 °F	1 Min	Yes	Shiny
Dentine 2	842 °F	6 Min	113 °F/Min	1670 °F	1 Min	Yes	Shiny
Glaze with glaze powder	1112 °F	2 Min	140 °F/Min	1670 °F	1 Min	---	Shiny

Alloys with a cte of 14.3 or higher may require a slow cool time of 2 – 3 minutes.

Non precious alloys may require a 3 minute cool time.

Follow alloy manufacturers instructions.

Horizontal furnace systems (Programat®) may require reducing the start temperatures to 752...