DES KROKODILS NÄCHSTE GENERATION
THE CROCODILE'S NEXT GENERATION

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With growing interest in low fusing glass Lucite ceramic systems, and having had success with Vision high fusing ceramic, Wohlwend AG was prompted to develop Vision Esthetic dual low fusing ceramic. Since its introduction at the 1999 IDS in Cologne, Germany, Vision Esthetic has established itself as an esthetically superior low fusing ceramic for both pressable and PFM restorations. Production processes are similar to those of Vision Classic, thus giving technicians the same advantages in shade accuracy and handling.

Dual Low Fusing Porcelain for Pressable and PFM restorations

- o *Compatible with PFM Alloys CTE range 15.9 - 16.7(25 - 500°C)*
- o *Kinder to natural - Vickers 470*
- o *Consistently accurate shade reproduction - V shades A1 - D4*
- o *Easy handling*
- o *No slumping during build up*
- o *Outstanding Esthetics*


The vital shade reproduction and natural fluorescence of Vision Esthetic corresponds to natural teeth. This replication of nature can be obtained by either using a basic two layer buildup or by incorporating additional dentine and incisal modifiers into the restoration. Vision Esthetic also offers a selection of 54 pressable ingots for build up and staining techniques, with opacities ranging from 25% to 85%.

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

Crowns:
Reduce the incisal edge by 2mm.
Labial / Buccal areas 1 - 1,5mm.
Distinct chamfer or 90° shoulder.
Palatinal minimum space of 1mm.
No sharp edges.

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Inlays / Onlays:

Inner surfaces should not Parallel, but open with an angle of 15° to the occlusal approx. 2mm deep.
No feather edges nor opposing bite at tooth / inlay junction.

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
Veneers / Laminats:

All prep. Margins as chamfers of 0.5 mm. Minimal thickness 0.5mm to the labial, anterior and up to 1mm to the incisal edge.

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Staining technique:

For inlays, onlays, veneers and rationally fabricated single crowns. After waxing to full contour, pressed with transparent ingots, the Objects are finished by using the appropriate stains.

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Layering technique:

Two ways to achieve high end aesthetics. Fully waxed up objects are pressed using the appropriate dentine ingot, Reducing the core as desired. This can also be done in wax. The core can now be characterized using stains and built up using dentins, Incisals and transpas.

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
Apply spacer in a thin layer 1mm above the prep. margin on crowns and veneers, to the margin on inlays and onlays.

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Stainig technique: Full conture wax up. No loss of dimensions during Pressing and firing processes.

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Build up technique: Either full wax up, reduced wax up Or waxed as a coping. Minimum thickness 0.8mm or 2/3 press material.

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Spruing: Wax wire 3-3,5mm Ø with a length between 4-7mm, should not taper toward the object Wax to the base in a slight angle taking care, that there is a distance of 3mm between objects.

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
Optional: Paper ring former and silicon ring former.

Remove strip from paper, form a ring to the appropriate mark. Slip onto base and adapt stabilizing ring. (These steps are not necessary when using siliconring). Pour special investment to just below the stabilizing ring. (Slightly over the mark in the silicon ring).

Carefully remove the stabilizing ring and fit the ring former, twisting it carefully. Investment must overflow through the opening. After bench set, remove paper / silicon ring, base and ring former, then clean the bottom surface with a plaster knife. This must be absolutely flat, so that the ring stands perfectly Upright in the pressing furnace. (IMPORTANT: Do not preheat the ingots!!!)

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Layering technique

(Optional also for staining)

Ingot	Opacity	Colour	Tendenz	Recommendation	Inlays / Onlays	Veneers	Crowns
A4 - D4	Approx. 40%	Accessed		A1 -D4	"X"	"X"	X
EA0 / EA00	Approx. 40%	Light, resp.	Vital preparation		"X"	"X"	X
EB0 / EB00	Approx. 40%	bleached shades			"X"	"X"	X
WA4 - WD4	Approx. 60%	Accessed	Covers slightly discoloured preparations	A1 - D4	"X"	"X"	X
WA0 / WA00	Approx. 60%	Light, resp.			"X"	"X"	X
WB0 / WB00	Approx. 60%	bleached shades			"X"	"X"	X
OP1	Approx. 85%	Light yellowish	Covers devital strongly stained Preparations	A1,A2,B1,B2,C1,C2,D2 A3, B3 A3,5, A4, B4 C3, C4, D3, D4	"X"	"X"	X
OP2	Approx. 85%	Redish yellow			"X"	"X"	X
OP3	Approx. 85%	Orange yellow			"X"	"X"	X
OP4	Approx. 85%	Kahki			"X"	"X"	X
OP5	Approx. 85%	White			"X"	"X"	X

"X" Layering technique, only on adequate preparation and clearance !

Staining technique

(Optional also for layering)

Ingot	Opacity	Colour	Tendenz	Recommendation	Inlays / Onlays	Veneers	Crowns
T1	Approx. 25%	Transparent			X	X	X
T2	Approx. 25-30%	Semi Transpa			X	X	X
CT1	Approx. 30-35%	Light yellowish	Mainly used for medium translucent teeth	A1,A2,B1,B2,C1,D2	X	X	X
CT2	Approx. 30-35%	Redish yellow		A3,A3,5,B3,D3	X	X	X
CT3	Approx. 30-35%	Orange yellow		A4,B4,D4	X	X	X
CT4	Approx. 30-35%	Greyish		C2,C3,C4	X	X	X
O2	Approx. 50%	Light whitish	For teeth with higher opacity levels		X	X	X
O2	Approx. 50%	Whitish			X	X	X

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Approximate value for calibrated press furnaces:

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

*Cut along the line,
using seperating disk.
Beware: Good suction!!!*

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
*Carefully sandblast objects by using glas beads. Start with 60 psi pressure, as soon as objects appear turn pressure down to 35 psi proceed blasting carefully.
DO NOT USE ALUMINIUM OXIDE*

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

Carefully separate the sprue from the objects with a sharp diamond disk. Check the fit and finish with fine diamonds burs and silicon rubber wheels. Do not use carbide burs or high speed handpieces.

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Prepare objects to desired contour (texture, contacts, etc) Sandblast with 50 μ aluminus oxide with low pressure and clean. Additions can be made, using overlay material. Customize the restoration according to situation, using Stains. This can be repeated several times. After desired effects are achieved, glaze once or twice as required. Firing is carried out on firing pillow.

Farbe / Colour	A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Opaker / Opaque	O	O	O	O	O	O	O	O	O	O	O	O	O	O	O	O
Opakdentin / Opaciuos body	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD
Dentin / Body	D	D	D	D	D	D	D	D	D	D	D	D	D	D	D	D
Schneide / Incisal	S	S	S	S	S	S	S	S	S	S	S	S	S	S	S	S
	57	58	58	60	60	57	58	59	60	58	58	59	60	57	58	59
Tranpamassen / Transpa materials	TR SCL CL OP RT BL WE NT GR															

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
Farbe / Colour	A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Schultermassen / Shoulder material	SP Neutral BA					SP Neutral BB				SP Neutral BC				SP Neutral BD		
Opaker zur Chromasteigerung Opaque to intensify chroma	OBA					O BB				O BC				OBD		
	10%	20%	30%	35%	40%	10%	20%	30%	40%	10%	20%	30%	40%	20%	30%	40%
Intensivopaker zum charakterisieren Intensiv opaques for characterization	I O WE OR BR V															
Dentin zur Chromasteigerung Body to intensify chroma	DBA					D BB				D BC				DBD		
	10%	20%	30%	35%	40%	10%	20%	30%	40%	10%	20%	30%	40%	20%	30%	40%
Intensivdentine zum charakterisieren Intensiv dentine for characterization	I D WE ST HG BR BL NEON															

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Thoroughly cleans the reduced, possibly individualized copings. Complete the reduced or missing portions, using the equivalent overlay materials. A correction cycle can be carried out if necessary. Prepare objects to desired contour (texture, contacts, etc) Finally glaze the restoration using stains to add small Nuances if desired.

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Notice minimum thickness of 0.8mm! Fire on firing pillow or thin metal pins, molars are only to be fired on firing pillow.

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Firing chart full-ceramic

Recommended Firing parameters!

Temperatures can vary due to different furnace performances

	Start temperature	Dry time	Heat rate	Final temperatur	Hold	Vacuum	Appearance
Staining / Internal staining technique							
Stain 1	842°F	4 Min.	131 °F/Min.	1436 °F	1 Min	---	
Stain 2	842°F	4 Min	131 °F/Min.	1436 °F	1 Min	---	
Glaze 1	842°F	4 Min	131 °F/Min.	1436 °F	1 Min	---	Slightly shiny
Glaze 2	842°F	Min	131 °F/Min.	1418 °F	1 Min		Shiny

Layering technique							
Dentine 1	842°F	6 Min.	131 °F/Min.	1436 °F	1 Min	Ja	Leicht Glänzend
Dentine 2	842°F	6 Min	131 °F/Min.	1418 °F	1 Min	Ja	Leicht Glänzend
Glaze without Glaze powder	842°F	4 Min	131 °F/Min.	1436 °F	1 Min	---	Shiny
Glaze with glaze powder	842°F	4 Min	131 °F/Min.	1400 °F	1 Min	---	Shiny

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Alloy Parameters

Compatible with alloys having a CTE of 15.9 – 16.7 (25° - 500°C).

Ideally at 16.2 (25° - 500°C)

If the CTE is Higher 16.4 (25° - 500°C)

- *Use a slow cooling cycle!*

If below 16.0 (25° - 500°C)

- *Remove restoration from lift as soon as furnace has fully opened!*

These parameters can vary between alloys!

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Metal processing prior to porcelain application

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
Process alloys according to manufacturers instructions (Delete "and avoid sharp edges"), Unless specified, the alloy should be sandblasted with min. 110µm AL²O² at 2 Bar.

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In most cases, especially on NP, the oxides are to be removed using min. 110µm aluminum oxide. (Refer to manufacturers manual!) A clean and grease-free surface must be present.

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*Apply opaque with a 75-80% semi-
covering layer.*

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
*Refer to firing chart.
Shiny appearance after firing.*

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*To achieve internal effects, characterizati-
ons can be made by using opaque-modi-
fiers "IO"*

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
To stabilize the shade in thin areas, apply a layer of opacious dentine.

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
*Incisal / Transpa (Optional.
To optimize light transmission)*

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After firing a shiny surface is desired.

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Firing chart PFM

Recommended Firing parameters
Temperatures may vary due to different furnace performances

	Start temperature	Dry time	Heat rate	Final temperatur	Hold	Vacuum	Appearance
Degassing	Follow alloy manufacturers manual						
Opaque 1*	842°F	6-8 Min.	131 °F/Min.	1472 °F	1 Min	Yes	Shiny
Opaque 2	842°F	6-8 Min	131 °F/Min.	1472 °F	1 Min	Yes	Shiny
Shoulder 1 / 2	842°F	4 Min	113 °F/Min.	1436 °F	1 Min	Yes	Shiny
Dentine 1	842°F	6 Min	113 °F/Min	1418 °F	1 Min	Yes	Shiny
Dentine 2	842°F	6 Min	113 °F/Min	1400 °F	1 Min	Yes	Shiny
Glaze with glaze powder	896°F	2 Min	113 °F/Min	1382 °F	1 Min	---	Shiny
Glaze without powder	896°F	2 Min	113 °F/Min	1418 °F	---	---	Shiny

***Final temperature for NP: 1562 °F**
Observe alloy CTE when considering a fast or slow cool time.

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To achieve opalescence in the incisal area, TR Opal can be mixed into all transpa materials. Even on a simple build up, natural light transmission and reflection can be achieved.

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
Mix glaze powder with liquid and apply in a thin layer. For individual characterisation use stains.

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Depending on desired level of translucency, apply TR Clear or Super Clear.

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Glaze appearance can be controlled by variation of hold time on final temperature.

WOHLWEND AG
DENTAL MANUFAKTUR

Internet www.wohlwend-ag.com . T 00423 373 42 43