WOHLWEND AG

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Introduction

The Vision-Classic is a high-fusing metal-ceramic and enables efficient and highly aesthetic work results. The outstanding features are the model lability and the extremely good stability of Vision-Classic.

The use of fluorescent and opalescent additional materials permits extremely individual and demanding results, optimally coordinated light dynamics and various possibilities, from simple dentin / incisal layering up to the masterly multilayering.

Compatibility with alloys in the CTE range from 13,8–15,2 $\mu\text{g}/\text{mk}$ (25–600°C)

- » High gold content or precious metal reduced
- » Non-precious metals
- » No green discoloration with silver-containing alloys
- » Highest colour accuracy with the V colours A1–D4
- » Easy handling
- » High stability
- » Little shrinkage
- » Excellent modelling properties

Classification

The Vision-Classic is a metal ceramic and a medical product of class IIa. The material is only for professional use in dental laboratories. Processing should only be carried out by trained personnel.

Important!

Observe alloy manufacturer information!

- » CTE of the alloy should be between 13.8-15.2 (25-600 °C).
- » Ideal value of the alloy is 14.1 to 14.2 (25-600 °C)
- » If the CTE value is higher than 14.5 (25-600 °C), long-term cooling is recommended.
- » If the CTE is lower than 14.1 (25-600 °C), the firing object is recommended to be removed immediately after the opening of the furnace chamber.
- » **In the case of non-precious metal alloys, only new metal may be used for casting the frameworks.**

The use of scrap metal can lead to a change of the CTE value and the risk of contamination of the alloy with foreign substances!

- » We recommend the use of Vision-Bonding.
- » Please follow the instructions of your alloy manufacturer!

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Modified Shade	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4	
Shoulder material	●	neutral + BA				neutral + BB				neutral + BC				neutral + BD			
Effect shoulder	●	ivory, straw, yellow-gold, brown, yellow-grey, brown-grey															
Opaque to intensify chroma	○	BA 10% 20% 30% 35% 40%				BB 10% 20% 30% 40%				BC 10% 20% 30% 40%				BD 20% 30% 40%			
Intensive opaques	○	white, yellow, orange, ocher, brown, violet, gingiva															
Body to intensify chroma	●	BA 10% 20% 30% 35% 40%				BB 10% 20% 30% 40%				BC 10% 20% 30% 40%				BD 20% 30% 40%			
Intensive dentine	●	ivory, straw, honeyyellow, brown, grey, blue, neon															
Intensive incisal	●	7, 8, 9, 10															
Mamelon	○	ivory, honey, sunset															
Effect masses	●	ivory, straw, honeyyellow, khaki															
Bleach	○	Opaque AB-0, Opaque-Dentin A0,B0, Dentine A0, B0, Incisal AB-0															
Cervical	●	straw, brown															
Gingiva	●	softpink, darkpink, softred, darkred															
Repair masses	○	Opaque, Dentine, Incisal, Transpa, Stain															

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





Metal processing prior to porcelain application

1. Prepare the alloys according to the manufacturer's instructions and round off sharp edges.
Unless otherwise indicated, the framework should be sandblasted with min. 110 µm aluminium oxide and 2 bar pressure. Make sure that the colour of the oxide is uniform. For non precarious alloys the oxide layer should again be radiated with min. 110 µm Aluminium oxide (follow the alloy manufacturer's instructions!). Clean surface and degrease with steamer or ultrasonic bath.
2. Apply Vision-Bonding Paste evenly and thinly. The following is fired according to bonder instructions.
3. The result after firing also replaces a wash firing.
4. Alternatively Vision-Classic Opaquer can be used for the wash firing. Apply the opaque thinly and evenly, see picture number 6.

	Start temperature Starttemperatur	Dry time Trocknen	Heat rate Anstieg	Final temperature Endtemperatur	Hold time Haltezeit	Vacuum Vakuum
Degassing / Oxydbrand	Follow alloy manufacturer manual / Legierungshersteller Angaben beachten					
Vision-Bonding	550°C / 1022°F	6 min.	80°C / 176°F min.	980°C / 1796°F	1 min.	Yes / Ja
Vision-Classic / Powder Opaquer 1	600°C / 1112°F	3 min	80°C / 176°F min.	970°C / 1778°F	1 min.	Yes / Ja
Vision-Classic / Powder Opaquer 2	600°C / 1112°F	3 min	80°C / 176°F min.	950°C / 1742°F	1 min.	Yes / Ja
Vis-Opak / Paste Opaquer 1	600°C / 1112°F	6 min	80°C / 176°F min.	970°C / 1778°F	1 min.	Yes / Ja
Vis-Opak / Paste Opaquer 2	600°C / 1112°F	4 min	80°C / 176°F min.	950°C / 1742°F	1 min.	Yes / Ja

1 st. Opaque:

Opaque powder & liquid

Mix the opaque to a creamy consistency with a brush or glass instrument than apply with 70% opacity. Firing temperature 970°C.

With the use of OL-2 Vision-Classic Opaque Liquid a similar viscosity as with Vis-Opaque Paste can be achieved.

Vis Opaque / Pastes Opaque

Apply with a brush or glass instrument. Firing temperature 970°C.

If necessary, the temperature can be reduced to the oxide firing temperature of the alloy. Holding time for 1 minute. Close tightly after use to prevent dehydration. Glossy surface after firing.

Second step: 2nd opaque.

1. apply the opaque even and covering.
2. slightly glossy surface after firing.

Options

To achieve effects characterized from the depth with the help of opaque modifiers „IO“

For example:

Place IO-Orange in the central fossa or in the cervical area with IO-Brown.

In order to reach a higher chroma from the depth, Chroma Opaque is applied to the equator or mixed according to colour chart.

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Schulterkombinationstabelle:

Third step: Ceramic shoulder / margin

Our shoulder masses in the basic colours A, B, C, D and Neutral will give you the possibility to create a shoulder for every colour situation. The mixing ratio can be taken from the table. Optionally, an even larger colour palette can be achieved with effect shoulder masses.

Shoulder/margin combination chart:

Vision – Classic Combination shoulder/margin chart																			
Shade	AO	BO	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4	
Shoulder/ Margin	Ivory + NT		Ivory + BA					Ivory + BB				Ivory + BC				Ivory + BD			
Mixing ratio in %	20		80	70	50	40	20	80	70	50	30	80	70	50	30	70	50	20	
	80		20	30	50	60	80	20	30	50	70	20	30	50	70	30	50	80	
Effect shoulder	ivory, straw, yellow-gold, brown, yellow-grey, brown-grey																		

	Start temperature Starttemperatur	Dry time Trocknen	Heat rate Anstieg	Final temperature Endtemperatur	Hold time Haltezeit	Vacuum Vakuum
Shoulder / Margin 1/2	450°C / 842°F	3 min.	45°C / 113°F min.	940°C / 1724°F	1 min.	Yes / Ja

Fourth step: layering

1. By limited space apply a thin layer of opaque-dentin (occlusal, cervical, etc).
2. Apply dentin (of the given shade) in full anatomical shape.
3. Create a small cut back in the incisal and approximal region. Build up the marginal ridges and mamelons with Transpa Neutral and mamelon masses (ivory, honey, sunset).
4. complete the tooth shape of the tooth with the corresponding incisal mass. A slight over dimensioning during the construction is ensured by the light Shrinkage during firing for a perfect result.
 - » The interdental space is filled with opaque-dentine.

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



Important:

If Vision-Classic ML modelling liquid is used, the starting temperature of the ceramic furnace must be observed. This should not exceed 450°C. The constituents of the modelling liquid need time to evaporate. If this time is not granted due to too high starting temperature and too low pre-drying time, the dry outer layer of the ceramic may cause it to an inner build up of steam, which can cause damage to the surface of your layered crown during the firing circle!

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





Fifth step: Correction

1. After the first Dentin firing.
2. If necessary, add shape corrections with Transpa materials.
3. Second dentin firing
4. After finishing (picture 4) the contact points, tooth shape and surface – macro and micro structures, proceed with staining and glaze firing. There are three glazing options available:
 - Glaze firing without glaze material and a subsequent mechanical glaze firing.
 - Polishing (rubber polishers, brushes, diamond polishing paste and so on).
 - Glaze firing with glaze material. Mix these with glaze liquid and apply thinly. If necessary, characterize with Visio-Classic Stains.
 - Glaze firing with Vision-Universal glaze and Stains.




Glaze Firing:



Layering Guide

- **Neck**
Straw / Brown / or
D.BA / BB / BC / BD
- Dentin
- TR-SCL
- ID.Honey
- Incisal / Enamel
- I-A0 / BO

- ID.Honey
- Dentin
- TR-SCL
- Incisal / Enamel
- I-A0 / BO



- **Neck**
Straw / Brown / or
D.BA / BB / BC / BD
- Dentin
- TR-NT
- MM 1 = ID.Ivory
- MM 2 = ID.Straw
- MM 3 = EF.Honey
- TR-SCL
- TR-BL
- II-10
- Incisal / Enamel
- II-7
- TR-NT

- TR-SCL
- MM 1 = ID.Ivory
- MM 2 = ID.Straw
- MM 3 = EF.Honey
- Dentin
- **Neck**
Straw / Brown / or
D.BA / BB / BC / BD
- TR-NT
- ID.Honey

- TR-NT
- **Incisal edge**
II-7
- Dentin
- **Light blocker**
OD.White & Dentin (1:2)
- TR-SCL
- **Mamelon**
MM 1 = ID.Ivory
MM 2 = ID.Straw
MM 3 = EF.Honey
- Incisal / Enamel
- II-10
- ID.Honey
- **Neck**
Straw / Brown / or
D.BA / BB / BC / BD

- Incisal
- Dentin
- Opak Dentin
- Opaker

Classic Firing chart

Guidelines! Deviations are possible due to different furnace performances.

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* With Vis-Opak paste, extend the drying time by 4 minutes to 6–7 pre drying time.
Observe alloy instructions! For horizontal furnace systems (e.g. Programat®)
Reduce start temperature to 400°C / 752°F
For REPAIRKIT please use separate firing chart!

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Wohlwend AG