LUKACryl SEP

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Release gel for 3D models

Discription

It is a medium viscosity gel that is applied to the surface of the 3D resin model with a brush and allowed to air dry. Despite its thickness, the gel spreads as soon as it is applied to the surface and forms a thin layer which, when dry, forms a thin, hard plastic layer that acts as a separating film between the model and the acrylic resin. It is effective after a single application. The thickness is 12 micrometers. A second layer can be applied for very rough surfaces.

Uses

LUKACryl Sep was developed to separate an orthodontic resin from a 3D resin printed model. Correct separation with other types of resins is not guaranteed, therefore a test should be carried out beforehand.

Advantages

- Simple application: a single coat is applied with a brush
- It forms a very thin release layer
- It does not contaminate the 3D resin or the acrylic
- It can be easily removed with water

Working instructions

Apply a few drops to the model and spread them with a brush. Rock for 30 seconds at low pressure in compressed air. For fast drying (30 sec) it is recommended to use compressed air with low pressure after application to avoid blowing. Otherwise, drying takes approx. 30 min. As it is not a lubricant, the acrylic resin layer must be removed from the 3D model with a knife tip after polymerization. The plastic film can be easily cleaned by washing with water and a brush.

It leaves no residue, does not affect the acrylic resin or the 3D model or change its properties and does not contain any toxic components.

Packaging

Bottle: 100 ml

Storage an expiry

Shelf life 3 years provided that the storage conditions are observed. Shelf life refers to the product in its original packaging.

Keep the product sealed in its original packaging in a dry place away from sunlight. After each use, the product should be properly resealed to prevent evaporation.

Warranty:

This information is based on the current state of our knowledge. Further information can be found on the safety data sheet of LUKACryl SEP.

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Lukadent GmbH