LUKAGingi/LUKAGingi Seperator

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INSTRUCTIONS FOR USE

LUKAGingi/LUKAGingi Separator

LABORATORY ADDITION SILICONE FOR THE CREATION OF REMOVABLE GINGIVAL MASK ON MODELS

DESCRIPTION / INDICATION

LUKAGingi: Medium consistency laboratory Vinylpolysiloxane for the reproduction of the alveolar mucosa. It is specific for works on models with implants. Being easy to trim, it allows a fast and accurate finishing touch. It can be scanned. Color: pink

LUKAGingi Seperator liquid separator to isolate LUKAGingi

WARNINGS

Before using LUKAGingi it is recommended to isolate the contact surfaces with LUKAGingi Seperator, in order to prevent adhesion with the materials of the same nature (example: addition silicone templates), or chemically different (condensation silicones or polyeters). Avoid contact with latex gloves.

DIRECTION FOR USE – VESTIGE GENGIVEST SEPARATOR

Vigorously shake the vial before use. Lay out LUKAGingi Seperator with the brush on the contact surfaces (dental impression or silicone masks), which must be clean and dry. Wait until Vestige gengivest separator is completely dry, facilitating this operation with an air jet.

PRECAUTIONS

Do not ingest, avoid contact with eyes

DIRECTIONS FOR USE – LUKAGingi

Insert the cartridge into the manual dispenser. The notch "V" on the basis of the cartridge should be facing down. Lock down the appropriate device of the dispenser. Remove the colored cap by rotating ¼ turn counterclockwise. By acting on the lever of the dispenser, dispense a small amount of material until the two components come out evenly. Delete the extruded material and mount the mixing tip in the reverse order to the removal of the cap. Proceed according to the chosen technique (direct or indirect).

CONSERVATION

Store in a dry place at a temperature of 5-27°C (41-80°F). Protect from heat. The materials are guaranteed until the expiration date when stored in their properly closed containers. The production lot and the expiration date are indicated on the packaging.

TECHNICAL DATA	LUKAGingi
Mixing ratio	1 :1
Mixing time (23°C/73°F)	Automixing
Working time (23°C/73°F)	2'
Setting time (23°C/73°F)	6'
Hardness - Shore A	72±2
Linear dimensional change	< 0,20%
Recovery from deformation	>99,8 %

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Lukadent GmbH