**LUKADENT<sup>®</sup>**

**LUKAVest 3D**

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## **Investment for 3D printing / light-curing burn-out modeling resins, crown and bridge technology as well as partial denture**

LUKAVest 3D is a phosphate-bonded, graphite-free, low-dust investment material especially suitable for all burn-out 3D printing resin and light-curing modeling resins. LUKAVest 3D is suitable for the entire crown and bridge technique as well as for partial denture. The rapid casting process (speed) only for waxed modelations. Conventional casting with oven preheating overnight should be used with light-cured resins. The corresponding expansion can be adjusted for each alloy type via the concentration of the corresponding Lukadent Liquids (see Concentration Recommendation).

Mixing Ratio:	100g powder: 24ml liquid 150g powder: 36ml liquid
Working time:	approx. 5 min
Recommended processing temperature:	19 °C

**Handling: Measure out the required amount of liquid and pour into the mixing bowl. Pour the investment into the same bowl. Mix manually for 20-30 seconds, then place in a vacuum mixing unit and mix for 1 minute.**

Bench set time:	approx. 20 min.
Furnace temperature:	max. 850 °C for rapid heating method
Compressive strength:	7-15 MPa

### **Preheat / burn out**

#### **Conventional (with holding levels) for light-curing resin**

- Fill casting ring with investment and allow to set. Allow the ring to cool to room temperature.
- Roughen the outer surfaces of the muffle and place in a cold furnace.
- Depending upon muffle size, allow to heat for 20-30 minutes at 250°C.
- Heat the furnace at a constant rate to final temperature.
- Allow the muffle to heat for a further 30-60 minutes.

#### **Rapid Heating for wax modelation**

- Fill the casting ring and allow to set for at least 20 minutes.
- Roughen the outer surfaces and place in a preheated furnace at final temperature.
- Allow the muffle to heat for 25-70 minutes at final temperature. Please consult the alloy manufacturers' instructions. (25 min size 1 ring; 30 min size 3 ring; 50 min size 6 ring; 70 min size 9 ring)

## Concentration table LUKAVest 3D

### Casting

Objects Alloy	Inlay 3 surface	Crowns & Bridges	Other telescope	Tap. Crowns 6°	partial denture
Precius > 70% Au	18 ml Liquid 18 ml Aqua dest.	18 ml Liquid 18 ml Aqua dest.	14 ml Liquid 22 ml Aqua dest.	14 ml Liquid 22 ml Aqua dest.	
Au reduced < 55% for ceramics	20 ml Liquid 16 ml Aqua dest.	20 ml Liquid 16 ml Aqua dest.	16 ml Liquid 20 ml Aqua dest.	16 ml Liquid 20 ml Aqua dest.	
Palladium based Alloys	24 ml Liquid 12 ml Aqua dest.	24 ml Liquid 12 ml Aqua dest.	20 ml Liquid 16 ml Aqua dest.	20 ml Liquid 16 ml Aqua dest.	
Non-precious Alloys		29 ml Liquid 7 ml Aqua dest.			18 ml Liquid 18 ml Aqua dest.

### Application examples

## Guarantee

All recommendations concerning application are based on our own experiences and tests and can be seen therefore only as approximate values. Our products are subject to continuous improvement. We reserve the right to change all parts and composition of this product. These instructions and recommendations are based upon current scientific standards and technology and can be considered correct to the best of our knowledge and experience. This version replaces all earlier versions.

Investment contains quartz and cristobalite. Avoid inhaling dust!

**Store the investment in a cool and dry place.**

---

Embedding is suitable for metal casting rings and for ringless systems. We recommend in combination with the investment, the use of specially matched products.

### **LUKACHrom C** **CrCo-Alloy for ceramics**

Composition in %  
Co 61,0 / Cr 25,0 / Mo 3,0 / W 8,0  
Vickers-Hardness 280 / CTE 13,9 / 14,0

### **LUKACHrom FH** **CrCo-Alloy for partial dentures**

Composition in %  
Co 62,5 / Cr 30,0 / Mo 5,0  
Vickers-Hardness 375


Casting Rings

Sprue Former

Casting Ring Liner

Manufactured of the highest quality pvc-material. Allows equal expansion in all directions. Produces precision castings. New and improved interlocking tapered ring & base would not allow accidentally separation. Transparency allows control of investment flow for a bubble-free environment. Easy-to-clean.

Lukadent GmbH  
Felsenbergweg 2  
D-71701 Schwieberdingen  
+49 (0) 7150 - 32955  
sales@Lukadent.de  
www.Lukadent.de

**VDDI**  
Dental Solutions.  
German Manufacturers.

Stand 12-2017