Felsenbergweg 2
D-71701 Schwieberdingen
+49 (0) 7150 – 32955
sales@Lukadent.de

LUKASil 85

Knetsilikon

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Anwendung

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Instruction for use

Indications for use

LUKASil 85 is an addition-curing 1:1 kneadable precision impression material with a final hardness of Shore A 85, its mixing ratio is very simple and save. LUKASil offers an extraordinary precise reproduction as a result of its detail graphic accuracy and a superior dimension stability due to a very low shrinkage, it is equipped with outstanding physical properties, high heat resistant, not sticky and physiologically uncritical.

LUKASil 85 is perfect suitable for check bites, silicone keys, bite registrations, models for denture repairs, blocking-out, teeth-insulating in the flask and all kinds of moulding. Working instructions

Working instructions

Take an equal part from both components (A white + B blue) and knead it intensively and in a timely manner approx. 30 sec. till a strip-free and homogeneous mixture is produced, recognizably by an evenly light blue colour. Thereupon the mass can be applied and after another approx. 8 min. it can be replaced.

Important working hints

For with-drawal and mixing gloves made from vinyl not latex should be used since latex extends the setting time or prevents it possibly. If mixing bare-handed the setting time will be poss. cut down

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Dosing:	1:1
Density:	1,60 g / cm ³
Colour code:	blau
Hardness:	Shore A 85
Mixing time:	ca. 30 sek.
Working time:	ca. 1,50 Minuten
Setting time:	ca. 8 Minuten
Lin. dim. change:	< 0,2%

Working room temperature
20°C – 23°C

Delivery form:
2 x 1,5 Kg tins white and blue
2 x 5,0 Kg buckets white and blue