**LUKADENT**

# LUKAGLAZE

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## Index

---

SAFETY.....	4
FEATURES AND TECHNICAL DATA.....	5
INSTALLATION INSTRUCTIONS .....	6
OPERATING INSTRUCTIONS .....	7-14
QUICK START INSTRUCTIONS.....	15
MAINTENANCE.....	16
TROUBLESHOOTING.....	16
ERROR CODES.....	17
PRODUCT SERVICE.....	18
WARRANTY.....	19
SUPPORT & CONTACT.....	20

# LUKA<sup>Glaze</sup>

## Glazing and Staining Furnace

---

The LukaGlaze is an economical, compact porcelain furnace for dentist clinic or laboratory air-fire staining and glazing. Final chairside esthetics are easily achieved with this portable, lightweight unit.

---

# SAFETY

---

## FOR PROFESSIONAL USE ONLY BY TRAINED DENTAL TECHNICIANS

Please read these operating instructions carefully before installing or operating this equipment.

- Use indoors only.
- Never operate the unit in close proximity to combustible materials or place materials on top of the unit, next to, or behind.
- The unit must be electrically grounded to a three wire electrical outlet or receptacle. The electrical service provided must be a dedicated line of the proper size according to local electrical codes.
- Unit must be placed in a position that allows the power cord to be easily disconnected from the wall or inlet socket.
- Do not attempt to maintain the unit until you have read and understand this operation manual.
- Disconnect the line cord before attempting to service the unit.
- Do not operate the unit controls with tongs or other tools.
- Do not use solvents or liquid cleaners on the control panel.
- Do not cover the top of the unit.
- The unit is not a burnout furnace or oven.
- If the unit is not operated in the manner as specified in this manual, the protection provided by the unit may be impaired.

## OSHA AND CALIFORNIA PROPOSITION 65: MUFFLE DUST EXPOSURE

In keeping with the policy of Lukadent GmbH to build safe products, comply with all National and State statutes and keep you, the valued customer informed; the services of a Certified Industrial Hygienist firm were employed to test and evaluate the lab operator's exposure to respirable refractory ceramic fiber (RCF) and crystobalite (a form of crystalline silica) present in the furnace muffle.

When it becomes necessary to replace the muffle, the person doing this work is recommended to wear a HEPA filter respirator and protective gloves as a precautionary matter.

Seal used muffle in a plastic bag and dispose of in accordance with local, state and Federal regulations.

Because this product and many similar products on the market today contain crystalline silica and ceramic fibers, it is necessary Proposition 65 that Lukadent includes the following statement:

"This product contains substance(s) to cause cancer."

## DESCRIPTION OF FURNACE AND INTENDED USE OF FURNACE:

The LukaGlaze is an economical, compact dental oven/furnace for dentist clinic or laboratory to air-fire of materials for dental prosthetics. The LukaGlaze is NOT intended for the burning out of binders, organic compounds and corrosive materials. The only intended use materials are dental stains and glazes.

# FEATURES & TECHNICAL DATA

---

## FEATURES:

- Small Footprint
- High Performance Muffle
- User-Friendly Interface with Dual Operating Modes
- Larger active heat zone for better temperature uniformity
- Full Temperature Range to 1100°C (2012°F)
- Smooth Door Movement
- Fast Heat Rates of up to 120°C/minute (216°F/minute)
- Energy Saver
- Easy temperature Adjustment
- Large Display with Easy to Read Parameters
- Limited warranty: 2 Year

## ELECTRICAL RATINGS:

Volts: 100/120VAC or 230VAVC, 50/60Hz

Watts: 600W @120V or 800W @230V

Watts to Maintain 1000°C (1832°F): 300W

## DIMENSIONS:

210mm (8.3") Width, 178mm (7.0") Depth, 222mm (8.8") Height

## WORK PLATFORM:

Honeycomb Tray: 55mm by 55mm (2.2" by 2.2"), 1.00mm Holes

## WEIGHT:

Furnace: 3.6Kg (8 Pounds) Shipping: 4.5Kg (10 Pounds)

## OPERATING CONDITIONS:

Temperature (Ambient): 5°C – 40°C (41°F – 104°F)

Relative Humidity: maximum 80%, non-condensing

## PERFORMANCE:

Temperature Accuracy: ±3°C from 400 to 1100C (±5.5°F from 752F to 2012°F)

Temperature Uniformity: ± 5°C (± 9°F) at steady state.

## RANGE OF FIRING PARAMETERS:

Programmable Value	Min		Max		Default	
	°C	°F	°C	°F	°C	°F
Low Temperature (Lo T), degrees	50	122	600	1112	500	932
Heat Rate (Rate), degrees/min	25	45	120	216	80	144
High Temperature (Hi T), degrees	601	1114	1100	2012	900	1652
Hold Time (Hold), minutes	0:00		3:00		1:00	

# INSTALLATIONS INSTRUCTIONS

## UNIT PACKAGE CONTENTS:

- One oven complete
- Power cord
- Owner & Operator's Manual (This document)
- Honeycomb Work Platform

## UNPACKING:


Save the box and other packaging material for future use in transporting the unit.

- Open the box.
- Remove packaging.
- Remove plastic bag.
- Remove all packaging material from inside the oven muffle and around the furnace.
- Unpack Honeycomb Tray and Power Cord

**ATTENTION:** Shipping damage should be reported to the carrier as soon as detected.

## INSTALLATION:

- Place Honeycomb Platform over rear ceramic rod then rotate down on front rod:
- Connect the furnace to a power circuit or receptacle with an over current protection of 15A or 20A.

The furnace should be located at least 15cm (6") away from walls, shelves and heat sensitive and combustible materials.

POWER SWITCH


POWER FUSE

POWER INLET

**BURN HAZARD IS PRESENT WHEN CHAMBER IS OPENED. USE CAUTION AFTER OPENING CHAMBER. INTERNAL AND EXTERIOR SURFACES MAY BE HOT.**

# OPERATING INSTRUCTIONS

## FUNCTION DESCRIPTIONS:

### POWER UP:


Turn on the power switch located on the back

At Power-up the LukaGlaze checks itself and performs an internal temperature calibration (About 60 seconds). It then goes into a Setup Routine where the user can adjust the Audio, Temperature Units, Temperature Offset and display Contrast.

**Audio:** Selection of alarm frequency and intensity  
(20 = high frequency, 0 = Mute)

Turn the Control Knob to change the value. Push Control Knob to store the new value. If the Control Knob is not turned pushing advances to the next setup parameter.


Only one setup parameter can be changed at a time. Turn off the oven and turn back on to adjust the next setup parameter.


# OPERATING INSTRUCTIONS

---

**°C/°F:** This parameter indicates which unit of temperature the control uses in operation. Turn the Control Knob CW or CCW to change units: Push Control Knob to store the new value.

**Tcal @ 960°C:** : Temperature offset Calibration of +/- 50°C is possible. Example: A value of +5°C increases the muffle temperature at 960°C by 5°C. A value of -5°C decreases the muffle temperature at 960°C by 5°C. Push Control Knob to store the new value


# OPERATING INSTRUCTIONS

---

The LukaGlaze finished the power-up and is now ready for operation. If the knob is NOT pushed then the control continues to end the Setup routine. The display shows "Done" for 2 or 3 seconds and then control proceeds to heat or cool to the Lo T and idle at this selected temperature. This is the Idle Mode when controlling at Lo T.

The first time the LukaGlaze is powered up the word "Test" will appear for 1 minute. In the "Idle Mode" the controller will heat up the muffle and control the temperature at the pre-programmed Lo T value (500°C). A heat-up test is performed for 1 minute.


# OPERATING INSTRUCTIONS

---


## **PROGRAMMING or REVIEWING PARAMETERS:**

In the Idle mode the 4 program parameters may be reviewed and changed. Parameters cannot be accessed once a program has been started.

A short push on Control Knob displays the Lo T value. Turn the Control Knob to change the Lo T value. Push the Control Knob to save the new Lo T value.

Two short pushes on Control Knob displays the Rate value. Turn the Control Knob to change the Rate value. Push the Control Knob to save the new Rate value

# OPERATING INSTRUCTIONS

---

Three short pushes on Control Knob displays the Hi T value: Turn the Control Knob to change the Hi T value . Push the Control Knob to save the new Hi T value

Four short pushes on Control Knob displays the Hold value . Turn the Control Knob to change the Hold value . Push the Control Knob to Save the new Hold value

# OPERATING INSTRUCTIONS

## SELECTION A PROGRAM / STARTING AND STOPPING:


A Long Push (>2 seconds) on Control Knob starts the Program at the current muffle temperature. The control will now heat to the preprogrammed Hi T at the programmed Rate value (see programming parameters). The approximate heating time will be displayed at the beginning of the second display line (4 minutes and 47 seconds).


Once the muffle temperature reaches the programmed Hi T the Hold timer (1:30) starts counting down.

The control will start beeping for 10 seconds before the end of the program. At the end of the programmed Hold time (1:30) the control will start beeping continuously to indicate the program is over. The program will terminate after a total of 3 minutes at the Hi T. A short push on the knob will turn off the sound but will not end the program.


Press and hold the knob (>2 seconds) to end the program and cool to Lo T.

# OPERATING INSTRUCTIONS

---

During the Rate or Hold turning the Control Knob will adjust the Hi T parameter; restart the heat rate and the hold time. If the operator has taken a quick look at the work and decided it is under-fired this is an easy way to continue the firing cycle with a modified Hi

If the knob is not pressed the program will end when the Hold time reaches 0:00 minutes. The display shows "End" and the control will cool the muffle to the Lo T as if the knob had been pressed

A pressing and hold will cancel the program at any time during a "Run". Pressing and holding the knob (>2 seconds) will cause the control to change between heating and cooling depending on the current muffle temperature and the programmed parameter values.

# OPERATING INSTRUCTIONS

---

## OPENING DOOR DURING IDLE:

When the control detects an open door condition the display will indicate Door. Pressing and holding the knob (>2 seconds) with the door open will prepare the unit to run a program. An alert will beep for 5 seconds if the door is not closed the control returns to Idle. If the door is closed during the 5 seconds the program starts.

## OPEN DOOR DURING RATE OR HOLD:


While the control heats the muffle during a program it is possible to open the muffle to check the work without canceling the program. However, the alert will sound for 5 seconds. After 5 seconds the program will be aborted and the muffle is turned off if the door has not been closed. If more time is needed to check the work, close the door after 3 or 4 beeps and open it again

# QUICK START INSTRUCTIONS

---

If the appliance is not used for 2 hours, a low temperature of 200 C is automatically set. By turning the control knob or opening the door, the oven heats up to the programmed Lo T again.

Turning the knob changes the programme. Press the knob to save the new programme.

The muffle is now heated or cooled with the programmed value to the displayed Lo-T of programme P-2.

# MAINTENANCE

---

## WARNING:

**This equipment contains dangerous voltages. Maintenance and repair work should only be performed by an authorized service technician of Lukadent GmbH. See Product Service on page 12 of this manual for contact information.**

To keep the equipment in good working order you should follow the guidelines below. Examine the equipment regularly for mechanical damage. Do not use solvents or liquid cleaners on the control panel.

## CLEANING:

Clean the unit as needed. Vacuum dust and dirt from the furnace rather than blow. This will minimize the amount of air born dust particles.

Clean exterior surfaces with a soft, damp cloth. Do not allow water to enter the control panel or around the knob.

**DO NOT IMMERSE IN WATER!**

# TROUBLESHOOTING

---

Problem	Diagnostic
No Power to the Unit	<ul style="list-style-type: none"> <li>- On/Off switch on the back of the oven is in the "0" off position.</li> <li>- Power cord is not fully plugged in to the back of the oven.</li> </ul>
Slow Temperature Rise	<ul style="list-style-type: none"> <li>- Check the programmed Rate value. Reprogram the Rate to a higher value.</li> <li>- Call Lukadent Customer Service and obtain Instructions on other diagnostic checks.</li> </ul>
Not holding Hi T	<ul style="list-style-type: none"> <li>- The maximum Hold time for the Hi T is 3 minutes. The control automatically cancels the program and cools to Lo T.</li> </ul>
Heating stops at Hi T when Door is open	<ul style="list-style-type: none"> <li>- The control automatically cancel the program is the door is open for more than 5 seconds max. The door can be opened repeatedly but each opening must be less than 5 seconds.</li> </ul>

# ERROR CODES

---

LCD Display	Diagnostic
Err 1 Over Temperature	The muffle temperature exceeded 1140°C Call Lukadent Customer Service
Err 2 Open TC	The thermocouple opened Replace thermocouple
Err 3 Open Muffle	The muffle heating wires opened Replace Muffle
Err 4 Wrong TC Connect	The thermocouple was incorrectly installed Reverse TC wires
Err 5 Tcal Error	A Temperature Calibration was incorrect Repeat Calibration
Err 6 EEPROM	Data could not be stored or read Call Lukadent Customer Service
Err 7 Check SD Card	The SD Card could not be read Call Lukadent Customer Service
Err 8 Unit Too Hot	The electronic components are too hot Cool the Unit
Err 9 Triac Short	The Unit cannot turn off current to the Muffle Call Lukadent Customer Service

# PRODUCT SERVICE

---

**SERVICE:**

Service of this equipment should be done only by qualified personnel. Failure to comply may result in personnel injury and/or property damage.

**METHODS OF PRODUCT SERVICE AVAILABLE:**

Telephone assistance available at the number listed below,  
Return the unit for servicing using the instructions below,  
Visit "www.lukadent.de" for additional information.

**BEFORE RETURNING THE UNIT:**

- Call prior to returning the oven. Tape a business card to the top of the oven to help identify problems with your oven.
- Do not ship the oven with the Honeycomb Working Tray inside the oven. Remove the Tray and pack in bubble pack outside the oven.
- Equipment damaged in shipment as a result of improper packing may not be paid by the carrier.

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**DISPOSING OF THE DEVICE:**

The device is an electronic device according to the "Act Governing the Sale, Return and Environmentally Sound Disposal of Electrical and Electronic Devices" (ElektroG).

It was identified in accordance with the existing law and provided with this symbol.

The device is not intended for private use. It is manufactured and delivered for commercial use and is to be disposed of by the end user according to the specifications of the Electrical and Electronic Equipment Act – ElektroG.

**CONFORMITY:**

Lukadent GmbH acknowledges that the following product:

Name: LukaGlaze™ was designed to conform to the basic requirements of the following EC guidelines:

Low Voltage Directive - Safety

EMC Directive - Electromagnetic Compatibility

Machinery Directive - Safety of Machinery

RoHS Directive - Regulation of Hazardous Substances

WEEE Directive - Waste from Electrical and Electronic Equipment reuse

# WARRANTY

---

**WARRANTY:** Except with respect to those components parts and uses which are hereinafter described, Lukadent warrants this furnace to be free from defects in material and workmanship for a period of two years from the date of sale. Lukadent liability under this warranty is limited solely to repairing or, at Lukadent option, replacing those products included within the warranty which are returned to Lukadent within the applicable warranty period (with shipping charges prepaid), and which are determined by Lukadent to be defective. This warranty shall not apply to any product which has been subject to misuse; negligence; or accident; or misapplied; or modified; or repaired by unauthorized persons; or improperly installed.

**INSPECTION:** Buyer shall inspect the product upon receipt. The buyer shall notify Lukadent in writing of any claims of defects in material and workmanship within thirty days after the buyer discovers or should have discovered the facts upon which such a claim is based. Failure of the buyer to give written notice of such a claim within this time period shall be deemed to be a waiver of such claim.

**DISCLAIMER:** The provisions here-in stated Lukadent sole obligation and exclude all other remedies or warranties, expressed or implied, including those related to MERCHANTABILITY and FITNESS FOR A PARTICULARPURPOSE.

**LIMITATION OF LIABILITY:** Under no circumstances shall Lukadent be liable to the buyer for any incidental, consequential or special damages, losses or expenses.

**LIMITATION OF ACTIONS:** The buyer must initiate any action with respect to claims under the warranty described in the first paragraph within one year after the cause of action has accrued.

Shipping Address:  
Lukadent GmbH.  
Felsenbergweg 2  
71701 Schwieberdingen  
+49 7150 32955  
www.Lukadent.de

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