MODE D'EMPLOI

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Dates techniques	Vestige putty soft fast	Vestige putty soft normal	Vestige putty ADmix	Vestige tray material	Vestige medium	Vestige light fast	Vestige light normal
ISO 4823/ADA spec. 19	Type 0 putty consistency	Type 0 putty consistency	Type 0 putty consistency	Type 1 heavy bodied consistency	Type 2 medium bodied consistency	Type 3 light bodied consistency	Type 3 light bodied consistency
Dosage Base: Catalyseur	1:1	1:1	5:1	5:1	1:1	1:1	1:1
Temps de mélange (23°C/73°F)	30"	30"	Automixing	Automixing	Automixing	Automixing	Automixing
Temps total de travail (23°C/73°F)	1'30"	2'	1'30"	1'30"	2'	1'30"	2'
Temps de prise en bouche (35°C/95°F)	2'30"	3'30"	3'	2'30"	3'30'	2'30"	3'30"
Temps de prise (35°C/95°F)	4'	5'30"	4'30"	4'	5'30"	4'	5'30"
Dureté Shore - A	64 ± 2	64 ± 2	64 ± 2	64 ± 2	50 ± 2	45 ± 2	45 ± 2
Variation dimensionnelle linéaire	< 0,2%	< 0,2%	< 0,2%	< 0,2%	< 0,2%	< 0,2%	< 0,2%
Mémoire élastique	> 99,5%	> 99,5%	> 99,5%	> 99,5%	> 99,5%	> 99,5%	> 99,5%

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Dati Tecnici	Vestige putty soft fast	Vestige putty soft normal	Vestige putty ADmix	Vestige tray material	Vestige medium	Vestige light fast	Vestige light normal
ISO 4823/ADA spec. 19	Type 0 putty consistency	Type 0 putty consistency	Type 0 putty consistency	Type 1 heavy bodied consistency	Type 2 medium bodied consistency	Type 3 light bodied consistency	Type 3 light bodied consistency
Rapporto di miscelazione	1:1	1:1	5:1	5:1	1:1	1:1	1:1
Tempo di miscelazione (23°C/73°F)	30"	30"	Automixing	Automixing	Automixing	Automixing	Automixing
Tempo totale di lavorazione (23°C/73°F)	1'30"	2'	1'30"	1'30"	2'	1'30"	2'
Permanenza nel cavo orale (35°C/95°F)	2'30"	3'30"	3'	2'30"	3'30'	2'30"	3'30"
Tempo di indurimento (35°C/95°F)	4'	5'30"	4'30"	4'	5'30"	4'	5'30"
Durezza Shore - A	64 ± 2	64 ± 2	64 ± 2	64 ± 2	50 ± 2	45 ± 2	45 ± 2
Variazione dimensionale lineare	< 0,2%	< 0,2%	< 0,2%	< 0,2%	< 0,2%	< 0,2%	< 0,2%
Recupero elastico	> 99,5%	> 99,5%	> 99,5%	> 99,5%	> 99,5%	> 99,5%	> 99,5%

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HIGH PRECISION POLYVINYLSILOXANE DENTAL IMPRESSION MATERIALS

DESCRIPTION / USE

Vestige putty soft fast / Vestige putty soft normal: putty consistency addition-curing silicone for the first impression, suitable for all techniques where hand-mixing material is required for filling the impression tray.

Vestige putty ADmix: putty consistency addition-curing silicone for impression taking with automatic mixing systems in 5:1 ratio.

Vestige tray material: heavy body consistency addition-curing silicone for impression taking with automatic mixing systems in 5:1 ratio.

Vestige medium: medium viscosity, hydro-compatible and thixotropic addition-curing silicone for high precision impression.

Vestige light fast / Vestige light normal: low viscosity, hydro-compatible addition-curing silicone for high precision impression.

Vestige products, combined according to the appropriate techniques, are suitable for impressions in fixed prosthesis, bridges, crowns, inlays, onlays, partial dentures, implants, position impressions, impressions for temporary and relining dentures.

WARNINGS AND PRECAUTIONS

Don't use in combination with polyether or condensation silicones. Avoid contact with latex gloves. Vinyl or nitril gloves are recommended. Carefully wash the hands if latex gloves have been worn, since interactions can be occurs (inhibition of the polymerization). Avoid contact with astringents containing sulphate. Carefully read the manufacturer's instructions. Do not use to patients with known or suspected sensitivity to silicone because it can cause allergic reactions. Do not swallow, avoid contact with eyes.

DIRECTIONS FOR USE

Vestige putty soft fast / Vestige putty soft normal – hand mixing

Take equal amounts of the two components with the proper scoop, being careful to manipulate only the material that will be used (to prevent contamination). Close the jars and start mixing until a homogeneous colour is obtained (around 30 secs). For a better retention to the impression tray, apply a thin layer of adhesive specific for A-silicone. Observe the manufacturer instruction for adhesive. Proceed on the base of the selected impression technique.

Vestige putty ADmix / Vestige tray material – 380 ml cartridge for 5:1 dynamic mixing systems

Preparation of the device:

Remove the cap-seal from the cartridge by grabbing the tab without excessive force. Insert the cartridge into the automatic mixing device following the instructions of the manufacturer of the device.

Apply the dynamic mixing tip and verify the correct insertion of the fitting hex drive shaft with that of the mixing tip. Then lock the device with the appropriate yellow plastic ring with a rotation of ¼ turn. If the operation is successful, at the start of the dispenser, the spiral inside the mixing tip must rotate. At first use, equalize the leakage of the two components by extruding a few centimeters of material and discard it or continue until the color of the material appears uniformly. After use, leave the dynamic mixing tip that will serve to cap until next use. Only use compatible dynamic mixing tips. To replace the cartridge, follow the instructions of the manufacturer of the automatic mixing device.

Vestige medium / Vestige light fast / Vestige light normal – 50 ml cartridge for 1:1 dispenser

Insert the cartridge into the manual dispenser 50 ml. The notch “V” on the basis of the cartridge should be facing down.

Lock down the appropriate device of the dispenser. Remove the cap by rotating ¼ turn counterclockwise.

By acting on the lever of the dispenser, dispense a small amount of material until the two components come out evenly. Delete the extruded material and mount the mixing tip in the reverse order to the removal of the cap. Insert the tip intra-oral for direct application on prepared teeth.

IMPRESSION TAKING

In case of the double impression technique (“two steps”), take the first impression with Vestige putty soft fast, Vestige putty ADmix or Vestige tray material. Once hardened, remove the tray from the mouth, rinse thoroughly and dry with spray air. With sharp instrument, remove undercut area from the impression as usual, then apply Vestige light fast (possibly also on the preparations) and proceed with taking the second impression. In the case of simultaneous impression technique (double mix or “one step”), apply Vestige light normal or Vestige medium directly on the preparations, simultaneously load the tray with Vestige putty soft normal, Vestige putty ADmix or Vestige tray material, thus introducing into the oral cavity and wait full setting.

DISINFECTION

After cleansing, the impression can be treated with specific disinfectants for elastomeric materials. Check and follow the instructions provided by the manufacturer.

PREPARATION OF MODELS / GALVANISATION

Impressions can be cast immediately after disinfection and they are not affected up to two weeks when stored properly. It is recommended to use model materials (plaster or resin), specific for dental use. The impression materials Vestige can be galvanized for electro-deposition and they are compatible with all CAD/CAM systems.

CONSERVATION

Store in a dry place at a temperature of 5-27°C (41-80°F). Protect from heat. The materials are guaranteed until the expiration date when stored in their properly closed containers. The production lot and the expiration date are indicated on the packaging.

Technical data	Vestige putty soft fast	Vestige putty soft normal	Vestige putty ADmix	Vestige tray material	Vestige medium	Vestige light fast	Vestige light normal
ISO 4823/ADA spec. 19	Type 0 putty consistency	Type 0 putty consistency	Type 0 putty consistency	Type 1 heavy bodied consistency	Type 2 medium bodied consistency	Type 3 light bodied consistency	Type 3 light bodied consistency
Mixing ratio	1:1	1:1	5:1	5:1	1:1	1:1	1:1
Mixing time (23°C/73°F)	30"	30"	Automixing	Automixing	Automixing	Automixing	Automixing
Total working time (23°C/73°F)	1'30"	2'	1'30"	1'30"	2'	1'30"	2'
Time in mouth (35°C/95°F)	2'30"	3'30"	3'	2'30"	3'30'	2'30"	3'30"
Setting time (35°C/95°F)	4'	5'30"	4'30"	4'	5'30"	4'	5'30"
Hardness Shore - A	64 ± 2	64 ± 2	64 ± 2	64 ± 2	50 ± 2	45 ± 2	45 ± 2
Linear dimensional change	< 0,2%	< 0,2%	< 0,2%	< 0,2%	< 0,2%	< 0,2%	< 0,2%
Recovery from deformation	> 99,5%	> 99,5%	> 99,5%	> 99,5%	> 99,5%	> 99,5%	> 99,5%

Higher temperatures reduce the times, lower temperatures increase them

Limitation of the responsibility: The indications related to the products TRAYART, both in oral form and written, are founded on the knowledge and experience acquired in dental field. The professional users are kept however to follow the indications for the use and to verify the fitness of the product according to the type of application. The professional user is assumed all the responsibilities for the consequential risks from a non correct use of the products TRAYART. A possible defectiveness of the products limits the responsibility of the manufacturer to the value of the products used by the professional.

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Technische Daten	Vestige putty soft fast	Vestige putty soft normal	Vestige putty ADmix	Vestige tray material	Vestige medium	Vestige light fast	Vestige light normal
ISO 4823/ADA 19	Type 0 putty consistency	Type 0 putty consistency	Type 0 putty consistency	Type 1 heavy bodied consistency	Type 2 medium bodied consistency	Type 3 light bodied consistency	Type 3 light bodied consistency
Dosierung Base : Katalysator	1:1	1:1	5:1	5:1	1:1	1:1	1:1
Anmischzeit (23°C/73°F)	30"	30"	Automixing	Automixing	Automixing	Automixing	Automixing
Gesamtverarbeitungszeit (23°C/73°F)	1'30"	2'	1'30"	1'30"	2'	1'30"	2'
Verweildauer im Mund (35°C/95°F)	2'30"	3'30"	3'	2'30"	3'30'	2'30"	3'30"
Abbindezeit (35°C/95°F)	4'	5'30"	4'30"	4'	5'30"	4'	5'30"
Härt Shore - A	64 ± 2	64 ± 2	64 ± 2	64 ± 2	50 ± 2	45 ± 2	45 ± 2
Lineare Dimensionsabweichung	< 0,2%	< 0,2%	< 0,2%	< 0,2%	< 0,2%	< 0,2%	< 0,2%
Ruckstellung nach Verformung	> 99,5%	> 99,5%	> 99,5%	> 99,5%	> 99,5%	> 99,5%	> 99,5%

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