MODE D’EMPLOI

VESTIGE IMPLA

IMPLANT IMPRESSION MATERIAL

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HIGH PRECISION POLYVINYLSILOXANE IMPLANT IMPRESSION MATERIAL

DESCRIPTION / USE

Vestige impla: implant impression material - heavy body consistency addition-curing silicone for impression taking with automatic mixing systems in 5:1 ratio. Indicated for indirect (closed tray) or direct (open tray) impression technique.

WARNINGS AND PRECAUTIONS

Don't use in combination with polyether or condensation silicones. Avoid contact with latex gloves. Vinyl or nitril gloves are recommended. Carefully wash the hands if latex gloves have been worn, since interactions can occur (inhibition of the polymerization). Avoid contact with astringents containing sulphate. Carefully read the manufacturer's instructions. Do not use to patients with known or suspected sensitivity to silicone because it can cause allergic reactions. Do not swallow, avoid contact with eyes.

DIRECTIONS FOR USE

Preparation of the device

Remove the cap-seal from the cartridge by grabbing the tab without excessive force. Insert the cartridge into the automatic mixing device following the instructions of the manufacturer of the device. Apply the dynamic mixing tip and verify the correct insertion of the fitting hex drive shaft with that of the mixing tip. Then lock the device with the appropriate yellow plastic ring with a rotation of ¼ turn. If the operation is successful, at the start of the dispenser, the spiral inside the mixing tip must rotate. At first use, equalize the leakage of the two components by extruding a few centimeters of material and discard it or continue until the color of the material appears uniformly. After use, leave the dynamic mixing tip that will serve to cap until next use. Only use compatible dynamic mixing tips. To replace the cartridge, follow the instructions of the manufacturer of the automatic mixing device.

DISINFECTION

After cleansing, the impression can be treated with specific disinfectants for elastomeric materials. Check and follow the instructions provided by the manufacturer.

PREPARATION OF MODELS / GALVANISATION

Impressions can be cast immediately after disinfection. It is recommended to use model materials (plaster or resin), specific for dental use. The impression materials Vestige can be galvanized for electro-deposition and they are compatible with all CAD/CAM systems.

CONSERVATION

Store in a dry place at a temperature of 5-27°C (41-80°F). Protect from heat. The materials are guaranteed until the expiration date when stored in their properly closed containers. The production lot and the expiration date are indicated on the packaging.

TECHNICAL DATA	VESTIGE IMPLA
ISO 4823 / ADA spec. 19	Type 1 heavy bodied consistency
Mixing ratio	5:1
Mixing time (23°C/73°F)	30'
Total working time (23°C/73°F)	1' 30"
Time in mouth (35°C/95°F)	2' 30"
Setting time (35°C/95°F)	4'
Hardness Shore - A	70 ± 2
Linear dimensional change	< 0,20%
Recovery from deformation	> 99,5%

Higher temperature reduce the times, lower temperatures increase them

VESTIGE BITE

BITE REGISTRATION VINYL POLYSILOXANE IMPRESSION MATERIAL

DESCRIPTION

Special formulated A-silicone for bite registration, thixotropic, fast and very hard after setting (95 Shore A).

INDICATION

All bite registration techniques, quadrant or full arch registration.

SETTING TIME

Total setting time: 1'30" from the start of the mix.

WARNING

Avoid contact with latex gloves. Vinyl or nitril gloves are recommended. Carefully wash the hands if latex gloves have been worn, since interactions can be occurs (inhibition of the polymerization).

PRECAUTIONS

Don't use toward patient with suspicious or verified sensibility to the silicone that can be cause of allergic reactions. Don't ingest, avoid contact with the eyes.

INSTRUCTIONS FOR USE

Insert Vestige bite in the manual dispenser for cartridges 1:1 from 50 ml. The nick to "V" on the base of the cartridge, it must be turned downward. Stop lowering the special device of the dispenser. Remove the colored cap rotating counterclockwise ¼ of turn in sense. Acting on the lever of the dispenser, extrude a small quantity of material up to that the two components escape in a uniform way. Eliminate the extruded material and insert the mixer performing the inverse operation to the removal of the cap. Progress applying Vestige bite on the occlusal plane and invite the patient to bite. Vestige bite is creamy and results imperceptible, so that not to jeopardize the naturalness in the occlusion. After the hardening, carefully remove the impression to avoid fractures.

DISINFECTION

After clearing, the impression can be treated with specific disinfectant for elastomeric impression materials. Please follow the manufacturer instructions.

PROCESS

Vestige bite can be finished with burs or sharpened instruments. It is dimensionally stable until two weeks. Vestige bite is scannable and compatible with the most used CAD/CAM systems.

CONSERVATION

Store Vestige bite in a dry place at temperature around 5-27°C (41-80 °F). The material is guaranteed until the expiration date when stored in properly closed packaging. The production lot number and the expiration date are indicated on the packaging.

TECHNICAL DATA	VESTIGE BITE
Mixing ratio	1:1
Mixing time (23°C/73°F)	Aurtomixing
Total working time (23°C/73°F)	30"
Time in mouth (35°C/95°F)	60"
Setting time (35°C/95°F)	1'30"
Hardness Shore - A	95
Linear dimensional change	< 0,05%
Detail reproduction	<20µm

Higher temperature reduce the times, lower temperatures increase them

Limitation of the responsibility. The indications related to the products TRAYART, both in oral form and written, are founded on the knowledge and experience acquired in dental field. The professional users are kept however to follow the indications for the use and to verify the fitness of the product according to the type of application. The professional user is assumed all the responsibilities for the consequential risks from a non correct use of the products TRAYART. A possible defectiveness of the products limits the responsibility of the manufacturer to the value of the products used by the professional.

FOR DENTAL USE ONLY

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VESTIGE IMPLA

IMPLANT IMPRESSION MATERIAL

VINYLPOLYSILOXAN ABFORMUNGSMATERIAL FÜR DENTAL ABFORMUNGEN ISO 4823/ADA spec.19: Type 1 heavy bodied consistency.

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