








TECHNISCHE DATEN · TECHNICAL SHEET · DONNÉES TECHNIQUE

		LUKASTONE IV SYNT-DIE Everest	LUKASTONE IV NATUR-DIE KZ	LUKASTONE IV RESIN-DIE Nanga	LUKASTONE IV EXTRA-HARD Makalu	LUKASTONE III SUPER-MODEL MontBlanc	LUKASTONE III ORTO KFO	LUKASTONE III MODEL Monviso	LUKASTONE II ARTI	LUKASTONE IV MODEL-BASE Mera	LUKASTONE II FLASK Anefo
ISO 6873 TYPE:CLASS		4	4	4	4	3	3	3	2:1	4	2:2
Druckfestigkeit 24h · Compressive strength · Résistance à la compression 24 h(MPa)		340	330	270	290	110	130	100	60	250	
Anmischverhältnis g / ml · Consistency · Consistance (ml)		20	20	20	22-24	30	27-28	28-30	30	24	50
Anrührdauer vakuum · Mechanical vacuum mixing · Mélange mécanique sous vide (sec)		30	30	30	30	30	30	30	-	30	-
Anrührdauer manuell · Manual mixing · Mélange manuel (sec)		60	60	60	60	60	60	60	30	60	60
Verarbeitungszeit · Working time · Temps de travail (min)		8	7	7	7	5	6	5	1.5	6	3
Entformungszeit · Final setting · Final temps de prise Vicat (min)		15	13	10	12	9	13	12	3	11	7
Expansion · Expansion · Expansion 2 h(%)		0.08	0.08	0.10	0.10	0.08	0.13	0.16	0.01	0.05	0.29
Druckfestigkeit nach 1h · Compressive strength · Résistance à la compression 1 h(MPa)		60	58	50	40	32	35	32	25	40	10

ANWENDUNGEN · APPLICATIONS SHEET · APPLICATIONS

Fester Zahnersatz Fixed prosthesis Prothèse fixe	CAD/CAM/Scannbar	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>								
	Implantate · Implants · Implants	Optional	Optional	<input checked="" type="checkbox"/>							
	Arbeitsmodell · Master cast · Maître modèle	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>							
	Zahnkranzgips · Dies · Dies	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>							
	Inlays-Onlays	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>							
	Brücken · Bridges · Bridges	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>						
Herausnehmbarer Zahnersatz Removable prosthesis Prothèse amovible	Gegenbissmodelle · Antagonists · Antagonistes (High Res.)	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>					
	Modellguss · Frameworks · Squelettique				<input checked="" type="checkbox"/>	Optional		Optional			
	Gegenbissmodelle · Antagonists · Antagonistes (High Res.)				<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Optional	<input checked="" type="checkbox"/>			
Kieferorthopädie · Orthodontics · Orthodontie	Kunststoffprothesen · Acrylic · Acrylique					<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>			Optional
					Optional		<input checked="" type="checkbox"/>				
					Optional		<input checked="" type="checkbox"/>				
Weitere · Other · Autres	Socket · Bases · Bases			Optional	Optional	Optional				<input checked="" type="checkbox"/>	
	Socket · Flask · Moufle							Optional			<input checked="" type="checkbox"/>
	Artikulator · Articulator · Articulateur								<input checked="" type="checkbox"/>		